



Shielding against X-Ray



marTECH Systems has expanded the product-range of shielding items.

X-Shield is a shielding material for X-Ray protection.

X-Shield is:

- usable everywhere
- simple to use
- economical
- available with a lead equivalent of 0,5 mm
- no special know how or tools required

Range of use:

- X-Ray units in hospitals or medical facilities,
in industrial or scientific installations



X-Shield - Installationsanleitung	X-Shield - installation guideline
<p>Untergrund: Der Untergrund muss vorbereitet sein wie zur Verklebung einer Dekortapete</p>	<p>Subsoil: The subsoil has to be prepared similar to the installation of a decorative wallcovering</p>
<p>Kleber: Verwenden Sie Bodenbelagskleber oder Dispersionskleber mit hoher Klebkraft, eine hohe Anfangshaftung wäre vorteilhaft</p>	<p>Adhesive: Use adhesive for flooring or other strong adhesive, high adhesion is needed</p>
<p>Verklebung: Kleber auf den Untergrund (im Allgemeinen Wände, Decke...) in geeigneter Menge auftragen. x-shield-Material in den Kleber einlegen und blasenfrei andrücken; die einzelnen Stücke werden an den Nähten jeweils 2 cm überlappend verklebt</p>	<p>Installation: Apply the adhesive to the subsoil (suitable quantity). Press the x-shield-material into the adhesive, no bubbles shall be left. At all joints overlap the material approx. 2 cm.</p>
<p>Erdung: Verwenden Sie selbstklebendes Kupferband, Artikel-Nr. 97406, und kleben Sie dies im Raum umlaufend auf die Oberfläche des x-shield-Materials. Verbinden Sie das Kupferband mit einer geeigneten Erdung</p>	<p>Earthing: Use self adhesive coppertape, item no. 97406, stick it on the surface of the x-shield-material, use the coppertape once around the room and connect it to a suitable earthing.</p>
<p>Achtung: Erdungsanschlüsse dürfen nur von autorisierten Fachkräften vorgenommen werden. Neben dieser Anleitung sind auch die Sicherheitshinweise für Handhabung, Verarbeitung und Entsorgung zu beachten</p>	<p>Warning: Earthing connections shall be done by authorised specialists only. Beside this guideline all valid standards and regulations have to be observed</p>

Shielding against X-Ray



X-Shield will be delivered in roles, compact and flexible



X-Shield is processed with usual tools and adhesives



The typical X-Shield installation consists of adhesive, X-Shield and a decorative wallcovering



Because of the overlapping adhesion will be arise a seam, which can be primed over by regular filler. The surface of X-Shield can be revised with a wallcovering or a coating.



Shielding against X-Ray

Details for Architects and Professional planner

You will get information for the necessary shielding from the manufacturer or the user of the X-Ray device.

X-Ray penetrates usual building materials, but alleviated. Hence, its possible that with building materials of high density and bigger material strength the use of an additional Shielding material is not necessary in any case. This case is often with floors and ceilings.

The Shielding effect / the required Shielding accomplishment is denoted a millimeter of thickness of lead, to the so-called „lead equivalent value“. The strengths required on this occasion are dependent on the used device. General details you find in German Institute for Standardization DIN 6812.

X-Shield is offered with a lead equivalent value of 0,5 mm.

If lead equivalent value with bigger value be required, so serveral layers of X-Shield can be used. Here, the single lead equivalent value were added.

As far as you wish our support with questions, texts or technical details, feel free to contact us:

marTECH Systems GmbH
Bertram-Schaefer-Strasse 11
35274 Kirchhain

Fon +49 6422-81-296

Fax +49 6422-81-228

info@marTECH-Systems.de

Shielding against X-Ray





X-Ray is used for many purposes in medicine and for technical matters.....

Too much X-Ray will cause health problems for humans

Humans have to be protected from higher doses of X-Ray

Manufacturers of X-Ray-equipment have to offer a radiation protection concept, which will be implemented in the construction of the building....





X-Shield is a simple way to fulfill the requirements of the radiation protection concept

The products name "X-shield" has been created according the the term "X-ray"





X-Shield ist:

- usable everywhere
- economical and simple to use
- available with a lead equivalent of 0,5 mm

- for different lead equivalent values a more layer construction can be used, the lead equivalents of the layer will be added



Range of use:

- hospitals
- medical facilities
- veterinarians
- industrial and scientific installations





... a decoration with a standard wallcovering or paint is possible





marTECH Systems shielding systems

- simple
- reliable
- economical

Prüfbericht Nr. LS-00605
Bestimmung des Bleigleichwertes nach DIN EN 61331-1

Auftraggeber: Marburg Technik, 35269 Kirchhain

Untersuchungsziel: Bestimmung Bleigleichwert bei 80 und 100 kV

Proben: Bleihaltiges Abschirmmaterial Produktname EMV X-Shield Nr. 97451 (0,35mm Pb-Nennwert), Nr. 97452 (0,50 mm Pb Nennwert)

Prüfdatum: 24.02.05

Prüfer: Dr. Eder

Prüfanordnung: Strahlengang gem. DIN EN 61331-1, Schmales Strahlenbündel, Messgerät PTW UNIDOS mit Luftkerma-Kammer M-23361. Es wurden jeweils 2 Proben gemessen und der Mittelwert gebildet. Die Abweichung der Proben untereinander war gering.

Messergebnisse:

Röhrensannung	Filter	Probe Art. Nr. 974541 Messwert (mm Pb)	Probe Art. Nr. 97452 Messwert (mm Pb)
80 kV	0,15 mm Cu	0,350	0,518
100 kV	0,25 mm Cu	0,351	0,514

Messunsicherheit +/- 3%

24.02.05



Dr.-Ing. Eder

Projectreport Marienhospital Altenessen, Essen

Shielding of a room for a Computer-Tomograph (CT) - Installation with X-Shield

Marienhospital Altenessen
Hospitalstrasse 24
45329 Essen, Germany

Thanks for a good cooperation to MEG Rhein-Ruhr, Mrs Berntßen



.... more than 100 years of tradition...





Marienhospital Essen

Projectreport

Orthopedic Clinic Munich-Harlaching

Shielding of a room for X-Ray use. X-Shield has been used for the installation

New building of the
Orthopedic Clinic Munich-Harlaching
Harlacher Strasse 51
81547 Munich

The installation has been done by: Weigert Corp., Weilheim





before installation



Orthopedic Clinic Munich-Harlaching



Preparation of the walls



Installation of X-Shield

Othopedic Clinic Munich-Harlaching



Installation of X-Shield



Orthopedic Clinic Munich-Harlaching



The surface of X-Shield is prepared with primer



A fibreglas wallcovering ist sticked on the surface of X-Shield

Orthopedic Clinic Munich-Harlaching



The fibreglas wallcovering is painted



The new room for X-Ray use is ready.

Projectreport

Clinik of the Friedrich-Schiller-University Jena

Shielding of a room for x-ray use by X-Shield

Thanks to Verbundgruppe Farbe & Heimtex GmbH, Jena for taking the pictures.

Clinik of the
Friedrich-Schiller-University
Bachstrasse 18
Jena, Germany



Clinik of the Friedrich-Schiller-University, Jena, Germany

Material and Tools



First steps in the installation



Cutting the material



Installation at a wall the the neighbour room



Apply the adhesive,
dispersionsadhesive or
adhesive for floor coverings



Installation of X-Shield





The installation is complete



A primer will make the installation of a decorative wallcovering more easy.



Clinic of the Friedrich-Schiller-University, Jena, Germany

X-Shield

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product name: X-Shield, item no. 97452

Manufacturer: marburg TECHNIC
A Division of Marburger Tapetenfabrik
J. B. Schaefer GmbH & Co. KG
P.O.Box 13 20
35269 Kirchhain / Germany

2. COMPOSITION/INFORMATION OF INGREDIENTS

Description: shielding of x-ray for walls, ceiling and floor,
lead containing foil

Hazardous components: None

Amount: Not applicable

EEC Symbol 231-100-4

CAR-No. 7439-92-1

Risk Phrases: According to 67/548/EWG apendix VI and
TRGS 505 solid lead containing material
does not offer risks for humans by breathing,
swallowing or skin contact

3. HAZARDS IDENTIFICATION

NOT classified as dangerous

Absorption of the materials by the body has to be avoided, especially breathing of
dusts / fractions as well as swallowing of material / dust or fractions.

Prevention-regulations have to be observed during processing.

4. FIRST AID MEASURES

Only necessary if material / dust / fractions or smoke are absorber by the body by
breathing or swallowing.

Eye: Washing for several minutes with plenty of water
Skin contact: Washing for several minutes with plenty of water
Inhalation: fresh air, in case of problems contact the doctor

5. FIRE FIGHTING MEASURES

Extinguishing media: CO2, Dry Chemical, Foam, Water
Special firefighting: in case of melted material do not use water
Fire & explosion hazards: Avoid inhalation of dust, steam and smoke
Autoplammability: Not available
Flash point: Not available

6. ACCIDENTAL RELEASE MEASURES

Spill and Leak procedures: Avoid inhalation of dust / steam / smoke,
use respirator mask,
avoid flow to gully / sewer

7. HANDLING AND STORAGE

Handling/storage: Keep dry,
no special requirements for the warehouse,
avoid exposition to acid

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Components: in case of use according to the guideline
no special protection required
Exposure limits: according to national regulations
Personal protection:
Glove Requirements: No special requirements,
washing of hands carefully after processing
is strongly requested
Eye protection: No special requirements
Clothing protection: No special requirements,
change of clothing after processing is needed,
wash the clothing carefully.

9. PHYSICAL AND CHEMICAL PRPOPERTIES:

Pure substance; Preparation or make:	Make
Physical form:	Solid
Colour:	silver / grey
Odur:	Non
Boiling point:	Not applicable
Melting/freezing point:	Not applicable
Sbility in water:	Not applicable
Autoflammability:	Not available
Flash point:	Not applicable
pH:	Not applicable

10. STABILITY AND REACTIVITY

Stability:	Stable – under normal application condition
------------	---

11. TOXICOLOGICAL INFORMATION

Effects of chronic esposure:	avoid absorption by the body – use of good industrial hygiene practices is recommended
Skin irritation:	None known
Hazard carcinogen :	None known

12. ECOLOGICAL INFORMATION

Persistence and degraddability :	Unknown
Bactiria toxicity:	Unknown
Fish toxicity:	Unknown
Bod:	Not applicable
Cod:	Not applicable

13. DISPOSAL CONSIDERATIONS

Waste Disposal methods:	Waste material has to be processed according to national regulations
Packing material warnings:	Not needed

14. TRANSPORT INFORMATION

IATA classification:	Not classified as dangerous
IMGD classification:	Not classified as dangerous
ADR/RID:	Not classified as dangerous
Hazardous substances Ordinance:	Not listed

15. REGULATORY INFORMATION

Not classified as dangerous

Directives 67/ 548 EEC and TRGS 505 have been considered when compiling this SDS; the information is provided for health and safety assessment by an industrial user. Reference should be made to any relevant local or national health, safety, and environmental legislation. This information does not constitute indication of suitability for specific uses.